

Veltek Associates, Inc. Healthcare Products for USP<797>

| VELTEK Item # | SIZE | DESCRIPTION | APPLICATIONS |
|--------------------|----------------|---|---|
| DECTR-08 | 24/16 oz | DECON-AHOL Trigger 16 oz -sterile | EPA Registered 70% IPA for all surfaces including, pass through items, compounding surfaces, hands/gloves, environmental surfaces, all ancillary equipment surfaces. Made with WFI and quadruple bagged. Wipe down of previously applied decontamination and cleaning agents. |
| DECWFI-SP-70 | 24/11 oz | 11oz aerosol mist -sterile | |
| DECWFI-SP-70-B | 24/11 oz | 11 oz aerosol inverta spray -sterile | |
| VEL12 12x12-S-3010 | 200 wipes | WIPEDOWN saturated 70% IPA wipe -sterile | |
| DQ200C-02 | 4/1 gal | DECON-QUAT200 C Concentrate -sterile | EPA Registered disinfectant, virucide, and fungicidal agent for use on all compounding and environmental (walls, floors, & ceiling) surfaces. |
| DQ200C-06-16Z-01 | 12/16 oz | Simple Mix RTU -sterile | |
| DQ200C-04-2Z | 4/1 gal | Simple Mix RTU -sterile | |
| SHC-02-5.25 | 4/1 gal | HYP0-CHLOR concentrate 5.25% -sterile | For use on inanimate, hard, non-porous environmental surfaces such as; walls, floors, ceilings, countertops, carts, stainless, etc. |
| SHC-16Z-0.52 | 12/16 oz | RTU .52% trigger spray -sterile | |
| SHC-16Z-0.25 | 12/16 oz | RTU .25% trigger spray -sterile | |
| VEL9-12x12-S-3025 | 100 wipes/case | Saturated 5.25% wipe -sterile | Ideal for deactivating hazardous cytotoxic compounds. |
| SPER-16Z-3% | 12/16 oz | STERI-PEROX 3% RTU trigger spray -sterile | For use on inanimate, hard, non-porous environmental surfaces such as; walls, floors, ceilings, countertops, carts, stainless, etc. |
| SPER-16Z-6% | 12/16 oz | 6% RTU trigger spray -sterile | |
| VEL10-12x12-3016 | 200 wipes/case | 6% saturated wipe -sterile | |
| DS200-02A | 4/1 gal | DECON-SPORE 200 Plus Concentrated -sterile | EPA Registered POAA sporicidal agent |
| DS200-07-16Z-01 | 12/16 oz | SimpleMIX RTU trigger spray -sterile | For use on all compounding and environmental (walls, floors, & ceiling) surfaces. Made with WFI and quadruple bagged. |
| DS200-04A | 4/1 gal | SimpleMIX RTU gallon pour -sterile | |
| DS200-03-13ZA | 12/13 oz | Concentrate unit dose | |
| DC-06-16Z-01 | 12/16 oz | DECON-CLEAN RTU SimpleMIX -sterile | Cleaning detergent to remove product and decontamination residues from all types of surfaces. Use prior to sporicidal application in rotation. |
| DC-04-1Z | 4/1 gal | RTU SimpleMIX -sterile | |
| VEL14 12x12-S-4003 | 200 wipes/case | Saturated wipe -sterile | |
| VEL13-9x9-S-3013 | 400 wipes/case | DRY-WIPE polycellulose blend 20x20 -sterile | Non-shedding drywipe for all decontamination & cleaning wipedowns |
| VEL13-9x9-NS-3024 | 1500/case | Polycellulose blend NS 300x5 | Hand drying – non sterile |
| STWA-02 | 4/1 gal | STERI-WATER USP Purified Water -sterile | For use with all sterile disinfectant concentrates |
| STWA-16Z | 12/16 oz | 16 oz Trigger spray -sterile | For rinse down of previously applied decontamination & cleaning agents |
| 1700-LC-1741 | Med - 25 | MITCOOL LAB COAT [disposable] -non-sterile | Heat –bonded polypropylene Lab Coat that is laminated with VAI's patented laminate binding material for ultimate barrier protection. Chemical resistant to a variety of agents |
| 1700-LC-1742 | Large - 25 | | |
| 1700-LC-1743 | X Large - 25 | | |
| 1700-LC-1744 | XX Large - 25 | | |
| 1700-PS-S-1792 | One size | Mitcool Protective Sleeves -sterile | Sterile sleeves for use when compounding sterile preparations. |
| PF-25M-2-S-2026 | One size | Face Mask w/ 2 Elastic Straps -sterile | Sterile face mask for use when compounding sterile preparations. |
| DH-09 | 12/32 oz | DECON-HAND -non-Sterile | Instant hand sanitizer before donning gloves. FDA approved. |
| DECWFI-BOT-01 | 12/32 oz | DECON-AHOL WFI -non-sterile 70% IPA | Glove disinfectant using 70% IPA with WFI. |
| DG-03-16Z | 12/16 oz | DECON-GLASS trigger RTU -sterile | Sterile glass and plexi-glass cleaner |
| SB-02 | 24/11 oz | STEEL-BRIGHT aerosol RTU -sterile | Sterile stainless steel cleaner |
| SBW-12x12-S -2315 | 100 1/pk | STEEL BRIGHT Saturated Wipes -sterile | |
| VEL13-123HD-S-3033 | 15 3/pk | WIPEDOWNHC- 1-2-3 Triple pack kits -sterile | Deactivation, neutralization, and decontamination of hazardous drug compounding surfaces. Easy 1-2-3 application of sterile saturated wipes. |
| | | | |

PERSONNEL GOWNING

VAI offers a complete line of sterile and non-sterile disposable garments to meet your needs. Our 1700 Mitcool blend is laminated with binding material for low particle shedding and provides excellent liquid chemical resistance. The material is also soft and breathable for optimal user comfort. The line includes:

- Coveralls
- Lab Coats
- Hoods
- Boots
- Sweat Bands
- Protective Sleeves
- Face Masks

CONSULTING AND TRAINING

Since 1981 VAI has been a leading global supplier and consultant of contamination control products serving the pharmaceutical, biotech, and healthcare facilities on cleanroom best practices. Our experienced team can help guide you and your facility on developing a sound cleaning and disinfectant program to help you meet the increasingly challenging regulations you face in your workplace today.

APPLICATION DEVICES

Our full line of product application equipment is designed for specific use with most all cleaning and decontamination products. All equipment and accessories are fully autoclavable and include:

- VAI's patented Core2Clean unit
- A variety of Mop and Bucket systems
- Sponges
- Covers

DISINFECTANTS, SPORICIDES, AND CLEANERS

VAI offers a full range of sterile and non-sterile cleaning agents, disinfectants, fungicides, virucides, sporicides, and sterilants including:

- 70% IPA
- Phenols
- Peroxyacetic Acid
- Hydrogen Peroxide
- Hypochlorite Solutions
- Quaternary Ammonium Compounds

All products are made with Water for Injection and are quadrupled bagged. All products come in a variety of packaging configurations and sizes including:

- RTU Trigger and Gallon Pour
- Aerosols (Spray and Mist)
- Gallon Concentrate
- Unit Dose Concentrates

CAPTURE CONTAMINATION

Sterile compounding procedures require clean facilities and air quality assessments. VAI offers a number of viable air monitoring solutions to ensure compliance with USP<797>. With our portable and built-in monitoring solutions, you can ensure that you are accurately and effectively exceeding the regulatory recommendations for sterile compounding in an aseptic environment.

ID CONTAMINATION

VAI LABS assist our customers by providing microbial testing services including but not limited to:

- Microbial Identification
- Antimicrobial Effectiveness Testing
- Time Contact Kill Studies
- Consulting

CLEANROOM WIPERS

VAI manufactures superior cleaning wipes for surface clean-ups, IPA wipedowns, Cleaning and Disinfecting of PEC's, stainless steel surfaces, glass, lexan, polyvinyles, and other critical surfaces required to be clean.

- 9"x9" sterile and non-sterile polycellulose dry wipes
- 12"x12" sterile premium polyester dry wipes
- Individually packaged sterile saturated 70% IPA (6"x6", 12"x12", and 18"x18")
- 12"x12" sterile saturated 70% IPA
- 3% and 6% sterile saturated H2O2 wipes
- .25 and .52% sterile saturated Hypochlorite wipes
- 5.25% sterile saturated Hypochlorite wipe for cytotoxic drug deactivation
- Individually packaged 12"x12" sterile saturated Steel-Brite wipes

